



DL2000 DNA Ladder

Product: DL2000 DNA Ladder, ready to use

Cat No: SG0530s

Package: 100 applications

Price:



1.5% Agarose, 5ul/lane
6cm lenth gel, 1xTBE, 15V/cm

Description:

Plasmid DNA was completely digested with appropriate restriction endonucleotides, ethanol precipitated and dissolved in storage buffer.

The DNA marker contains the following fragments: 100, 250, 500, 750, 1000,2000bp.

Quantity per Package:

100 applications.

Suggest Loading Amount:

5ul/lane.

Storage & Loading Buffer:

10mM Tris-HCl (pH7.6), 10mM EDTA, 0.01% xylene cyanol FF, 0.01% Bromophenol Blue, 15% glycerol.

Storage:

4°C

Note:

If the width of Agarose gel lane is more than 6mm, please load more ladder DNA onto the gel.

Do not heat before use.

Suggest gel concentration: 1.0-2.0%.



上海闪晶分子生物科技有限公司

地址：上海市闵行区北桥镇吴河路328号A栋2楼

邮编：201109

联系：市场部

电话：54460832 800-988-1995

E-mail:master@shinegene.org.cn

网址：www.shinegene.org.cn

